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The reduced horizontal crest volume is a very frequent clinical problem. Atrophy causes problems both from a biological and a surgical point of view, that we must evaluate in clinical choices. The thin atrophic ridges limit the possibility of inserting dental implants without first performing reconstructive procedures that increase surgical invasiveness and post-operative morbidity, and the standard positioning of screw implants in thin ridges often results in a vestibular dehiscence that represents the main anatomical cause of peri-implantitis and the related implant failures. For this reason, thin crests require bone reconstruction with grafts or GBR. Alternatively, the patient with horizontal bone atrophies can be treated with special wedge implants for thin ridges. In a clinical case, we have applied the principles of minimally invasive implant surgery which can be obtained using cuneiform implants inserted with press-fit technique. These implants use a site prepared with piezosurgery and insertion with magnetic thrusts.

BIBLIOGRAFIA

1. Amner MH. The time sequence of tissue regeneration in human extraction wounds. *Oral Surg Oral Med Oral Pathol* 1969; 27: 309-318.
2. Johnson K. A study of the dimensional changes occurring in the maxilla following tooth extraction. *Aust Dent J* 1969; 14: 241-254.
3. Schropp L, Wenzel A, Kostopoulos L, Karring T. Bone healing and soft tissue contour changes following single-tooth extraction: a clinical and radiographic 12-month prospective study. *Int J Periodontics Restorative Dent* 2003; 23: 313-323.
4. Araújo MG, Lindhe J. Dimensional ridge alterations following tooth extraction. An experimental study in the dog. *J Clin Periodontol* 2005; 32: 212-218.
5. Spinato S, Rebaudi A, Bernardello F, et al. Piezosurgical treatment of crestal bone: quantitative comparison of post-extractive socket outcomes with those of traditional treatment. *Clin Oral Implants Res*. 2016 Mar;27(3):361-6.
6. Gupta A, Tiwari B, Goel H, Shekhawat H. Residual ridge resorption: a review. *Indian J Dent Sci* 2010; 2:7-11.
7. Tallgren A. The continuing reduction of the residual alveolar ridges in complete denture wearers: a mixed-longitudinal study covering 25 years. *J Prosthetic Dent* 2003; 89: 427-435.
8. Milinkovic I, Cordaro L. Are there specific indications for the different alveolar bone augmentation procedures for implant placement? A systematic review. *Int J Oral Maxillofac Surg* 2014; 43: 606-612.
9. Esposito M, Grusovin MG, Felice P, et al. The efficacy of horizontal and vertical bone augmentation procedures for dental implants - a Cochrane systematic review. *Eur J Oral Implantol* 2009; 2: 167-184.
10. Vercellotti T. Piezoelectric surgery in implantology: a case report - a new piezo-electric ridge expansion technique. *Int J Periodontics Restorative Dent* 2000; 20: 358-365.
11. Azi ML, Aprato A, Santi I, et al. Autologous bone graft in the treatment of post-traumatic bone defects: a systematic review and meta-analysis. *BMC Musculoskelet Disord* 2016 Nov 9;17(1):465.
12. Ronda M, Rebaudi A, Torelli L, Stacchi C. Expanded vs. dense polytetrafluoroethylene membranes in vertical ridge augmentation around dental implants: a prospective randomized controlled clinical trial. *Clin Oral Implants Res* 2014 Jul;25(7):859-66.
13. Linkow LJ, Winkler S, Shulman M, et al. A new look at the blade implant. *J Oral Implantol* 2016 Aug;42(4):373-80.
14. Linkow LJ. The blade-vent: a new dimension in endosseous implantology. *Dent Concepts* 1968; 11: 3-12.
15. Linkow LJ, Cherchéve R. Theories and techniques of Oral Implantology vol. 1, St. Louis, MO: Mosby, 1970.
16. Marcello-Machado RM, Faot F, et al. Mini-implants and narrow diameter implants as mandibular overdenture retainers: A systematic review and meta-analysis of clinical and radiographic outcomes. *J Oral Rehabil* 2018 Feb;45(2):161-183.
17. Vercellotti T, Nevins ML, Kim DM, et al. Osseous response following resective therapy with piezosurgery. *Int J Periodontics Restorative Dent* 2005; 25: 543-549.
18. Preti G, Martinasso G, Peirone B, et al. Cytokines and growth factors involved in the osseointegration of oral titanium implants positioned using piezoelectric bone surgery versus a drill technique: a pilot study in minipigs. *J Periodontol* 2007; 78: 716-722.
19. Stacchi C, Vercellotti T, Torelli L, et al. Changes in implant stability using different site preparation techniques: twist drills versus Piezosurgery. A single-blinded, randomized, controlled clinical trial. *Clin Implant Dent Relat Res* 2013; 15: 188-197.
20. Vercellotti T, Stacchi C, Russo C, et al. Ultrasonic implant site preparation using Piezosurgery: a multicenter case series study analyzing 3,579 implants with a 1- to 3-year follow-up. *Int J Periodontics Restorative Dent* 2014; 34: 11-18.
21. Szmukler-Moncler S, Salama H, Reingewirtz Y, Dubruille JH. Timing of loading and effect of micromotion on bone-dental implant interface: review of experimental literature. *J Biomed Mater Res*. 1998
22. Reilly DT, Burstein AH. The elastic and ultimate properties of compact bone tissue. *J Biomech* 1975; 8: 393-405.
23. Guo EX. Mechanical properties of cortical bone and cancellous bone tissue. In: *Bone Mechanics Handbook*, Cowin SC, ed. Boca Raton, FL: CRC Press, 2001, pp. 10/1-10/23.
24. Abuhussein H, Pagni G, Rebaudi A, Wang HL. The effect of thread pattern upon implant osseointegration. *Clin Oral Implants Res* 2010 Feb;21(2):129-36.

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© RIPRODUZIONE RISERVATA